

A close-up photograph of a laser dilatometer probe. The probe is a clear, cylindrical glass tube with a smaller inner tube. It is mounted on a polished, metallic base. The background is a blurred, light blue surface. The entire image has a blue color cast.

## Laser Dilatometer L75

**LINSEIS**




Normally available dilatometers of the series L75 which use the push rod principle of measurement are successfully available in the market since 1957. As a further total new development Linseis is now offering a Laser Dilatometer built after the Michelson-Interferometer principle.


The original push rod dilatometer is based on a development made by Dr. Max Linseis in the early 1950<sup>th</sup>. In those days the principle used was the same as today, using a measuring tube in which the sample is placed. The expansion signal of the sample is transmitted through a push rod onto a LVDT sensor (linear variable differential transformer). This sensor has theoretically an indefinite resolution. The push rod, which transmits expansion Together with other system components and the following evaluation electronics a maximum resolution of 10nm can be obtained.

The push rod, which transmits expansion signal from the sample to the sensor, always covers a temperature gradient of up to 1600°C from end to end, sometimes even up to 2400°C. It is understandable that these high temperatures differences on these short dimensions makes it very important for good reproducible results, to insure very good thermal conditioning of the push rod / measuring system. If this expansion transmitting part will see a different thermal history from run to run, the reproducibility can not be best. In order to get best results, identical conditions from run to run are a must.

For example the room temperature, the start temperature, the heat up speed and the used sample atmosphere etc. have all to be identical.



A further increase of resolution, as well of the absolute accuracy is now also possible through the new development of the LINSEIS Laser Dilatometer of the Pico-series. As the name indicates already the resolution goes up to Picometers (0,3nm = 300 Picometer). That means resolutions can be obtained which are up to a factor 33,33 higher as the resolution that were possible up to date. On top the principle of interference measurement give the possibility for much higher accuracy• s, especially as some special computer calibrations are used. Up to now absolute accuracy• s of 1% were normal, with best accuracy• s up to 100nm. The new method allows accuracy• s up to 30nm.



<b>Range</b>	20 mm
<b>Resolution</b>	± 1 nm (±0.3 nm typical)
<b>Accuracy</b>	± 50 nm
<b>Reproducibility</b>	± 10 nm
<b>Wave length</b>	632.8 nm
<b>Ambient temperature</b>	10 to 30°C
<b>Sample holder</b>	Fused silica
<b>Push rod</b>	Fused silica, evacuated
<b>Heating</b>	Inductive
<b>Frequency range generator</b>	150 to 400 kHz
<b>Temperature range</b>	RT to 1000°C
<b>Thermocouple</b>	Type-K (NiCr/Ni), welded to sample
<b>Heating/cooling rate</b>	0.1 to 99.9 K/s / 0.1 to 99.9 K/min
<b>Dwell time</b>	0 to 3000 s / 0 to 3000 min
<b>Acquisition rate</b>	max. 1000 readings/s
<b>Power requirements</b>	230 VAC / 16 A / 50 to 60 Hz
<b>Power consumption generator</b>	3.5 kVA
<b>Cooling water</b>	2 l/min / 5 bar
<b>Sample/cooling gas</b>	Inert gas, 10 l/min max. at 1 bar max.



<b>LINSEIS GmbH</b>	<b>LINSEIS Inc.</b>
Vielitzerstr. 43	109 North Gold Drive
95100 Selb	Robbinsville, NJ 08691
Germany	USA
Tel.: (+49) 9287-880 - 0	Tel.: +01 (609) 223 2070
Fax: (+49) 9287-70488	Fax: +01 (609) 223 2074
E-mail: <a href="mailto:info@linseis.de">info@linseis.de</a>	E-mail: <a href="mailto:info@linseis.com">info@linseis.com</a>

**Products:** DIL, TG, STA, DSC, HDSC, DTA, TMA, MS/FTIR, In-Situ EGA, Laser Flash, Seebeck Effect  
**Services:** Service Lab, Calibration Service

[www.linseis.com](http://www.linseis.com)